
Biopsie-Ventil ohne Spülkanal SU

DE

EN

FR

IT

ES

NL

CZ

PT

SE

NO

DK

RU

CN

KR

TU

PL

SA

EL

HR

RO

DE	Biopsie-Ventil ohne Spülkanal SU	Gebrauchsanweisung	3
EN	Biopsy valve without flushing channel SU	Instruction Manual	8
FR	Valve de biopsie sans canal de purge su	Manuel d'utilisation	13
IT	Valvola per biopsia senza canale di spurgo su	Manuale di istruzioni	18
ES	Válvula de biopsia sin canal de purga su	Manual de instrucciones	23
NL	Biopsieklep zonder spoelkanaal SU	Handleiding	28
CZ	Bioptický ventil bez proplachovacího kanálu SU	Návod k použití	33
PT	Válvula de biópsia sem canal de irrigação SU	Manual de Instruções	38
SE	Biopsiventil utan spolningskanal SU	Instruktionsbok	43
NO	Biopsiventil uten spylekanal SU	Instruksjonshåndbok	48
DK	Biopsiventil uden vandingskanal SU	Betjeningsvejledning	53
RU	Биопсийный клапан без ирригационного канала SU	Руководство по эксплуатации	58
CN	无灌注通道的活检阀SU	说明书	63
KR	플러싱 채널 SU가 없는 생검 밸브	사용 설명서	67
TU	Yıkama kanalı SU olmayan biyopsi valfi	Kullanım talimatları	72
PL	Zawór bioper bez kanału płuczającego SU	Instrukcje użytkowania	77
SA	SU فيظنت فانق رودب ءعزخا امامص	مادختس الل تاميلعت	82
EL	Βαλβίδα βιοψίας χωρίς κανάλι άρδευσης SU	Οδηγίες χρήσης	86
HR	Ventil za biopsiju bez ispiranja kanala SU	Uputa za uporabu	91
RO	Valvă de biopsie fără canal de irigare SU	Instrucțiuni de utilizare	96

Symbole

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Symbols

EN

	Medical Device		Latex Free
	Catalog Number		Non-Sterile
	Batch		
	Warning! Read the Instructions for use before use		
	Usable until		
	For single use only. Do not reuse		
	Do not re-sterilize!		
	Keep dry		
	Protect from sunlight		
	Manufacturer		
	Date Of Manufacture		
	Do not use if packaging is damaged		

Copyright

The copyright to these instructions for use remains with the company ENDO-FLEX GmbH. All rights, including those of translation, are reserved. No part of this information may be reproduced in any form (print, photocopy, microfilm or any other process), in whole or in part, without the written permission of the company ENDO-FLEX GmbH or processed, reproduced, translated or distributed using electronic systems. Violations can result in criminal penalties.

Limitation Of Liability

The does ENDO-FLEX GmbH not assume any liability for personal injury, damage to property, damage to the product or consequential damage caused by non-compliance with these instructions for use, improper use of the product or any other action by unqualified and not ENDO-FLEX GmbH trained personnel. It is also forbidden to carry out modifications or technical modifications to the product by yourself.

Purpose of the Instructions for use

These Instructions for use are intended for use by the medical personnel working with the product. It describes the structure, function and operation of the product. Read these Instructions for use carefully and make sure that you understand them and that they are accessible to healthcare professionals do it!

Follow the instructions to ensure safe working. Also, follow the instructions in other devices used in combination with this product.

Do not follow instructions from unqualified third parties. The use of the product may only be carried out by qualified and, if necessary, on request ENDO-FLEX GmbH by trained medical personnel. The use of this product by healthcare professionals is always carried out taking into account the individual benefits of the patient at their risk.

Classification

The biopsy valve is designated as a Class I product according to Medical Device Regulation (MDR) – (2017/745) on medical devices.

Reporting Obligation

The indicates ENDO-FLEX GmbH that the user is obliged to report to the competent authority of his country and at the same time all adverse events occurring during the use of this product and concerning this ENDO-FLEX GmbH product.


Language

These and all other instructions for use are available as a PDF download at www.endo-flex.de .

1. Product Description

The biopsy valves offered by ENDO-FLEX GmbH consist of an elastic and non-slip plastic, which facilitates handling, especially with protective gloves.

Color differences (red and blue) almost eliminate confusion with compatible endoscope types. The removable cap also allows the use of instruments with larger external diameters.

1 Biopsy Valve

This Instructions for use is valid for products listed below:

Biopsy valves without flushing line

- ZB0100
- ZB0120

1.1. Contents of the packaging

- 1 Carton
- 100 Biopsy valves SU (Single use) individually non-sterile packed
- 1 Instructions For Use

1.2. Intended Use

The biopsy valve covers the biopsy/intake channel inlet of a gastrointestinal endoscope. The biopsy valve enables the introduction and exchange of endoscopic instruments, supports the maintenance of insufflation and minimizes leakage during gastrointestinal endoscopy from the biopsy channel inlet.

1.2.1 Duration Of Use

The application duration is specified for less than 60 minutes.

1.2.2 Lifetime

The life of the product is 3 years.

1.2.3 Indication

- Endoscopic Examinations

1.2.4 Contraindications

- No contraindications are known for this product.
- Possible injuries associated with endoscopic examinations may include: Perforations, bleeding, infections such as acute pancreatitis.

1.2.5 Potential Complications/Risks

- Product-related complications and/or risks are not known.
- Sedation during endoscopic examination increases the risk of hypoxemia, hypercapnia, hypotension, arrhythmias and aspiration due to reduced protective reflexes. Hypoxemia also occurs without sedation in endoscopic examinations due to the intake of the endoscope.
- Possible injuries associated with endoscopic examinations may include: Perforations, bleeding, infections such as acute pancreatitis.

1.2.6 Target group of users

Trained specialists and specialists.

1.3. Patient Population

The patient or patient target group is derived from the indication of the responsible physician, who treats the patient in a diagnostic or therapeutic manner within the scope of an endoscopy (the leading procedure in itself), according to the intended use of the medical device. Restrictions on patient populations and/or target groups are not known.

1.3.1 Use of the product in the case of minors

The indication for the use of the product in minors must be limited by the user on the basis of the respective individual physiological and anatomical conditions.

1.3.2 Use of the product in pregnant or breastfeeding women

The indication for the use of the product in pregnant or breastfeeding women must be limited by the user on the basis of the respective individual physiological and anatomical conditions.

1.4. Sterility

These disposable medical devices are supplied non-sterile. Treatment and sterilization is no longer possible after use and is prohibited! The product may only be used once and must be disposed of after use according to national and internal guidelines.

2. Safety Information

The safety and warning information used in these Instructions for use is based on ANSI Z535 standards.

DANGER	Warns of an imminent danger that leads to death or serious injury.
WARNING	Warns of a potentially dangerous situation that could lead to death or serious injury.
CAUTION	Warns of a potentially dangerous situation that could lead to moderate or minor personal injury.
NOTE	Warns of a potentially dangerous situation that could lead to property damage and environmental damage.
SAFETY INSTRUCTIONS	Safety instructions (or equivalent signs) identify specific safety-related instructions or procedures.

2.1. General

SAFETY INSTRUCTIONS

Take appropriate measures to protect the patient and yourself from contamination and injury. Follow national and internal hygiene and safety guidelines. If in doubt, contact your local representative.

WARNING

Do not use defective instruments! If any defects occur, discard the instrument and replace it with a new one.

NOTE

Never use the product outside the recommended technical specifications.

WARNING

Never manipulate the structural conditions of the instrument; in the event of a malfunction, discontinue use.

3. Application

3.1. Work Preparation & Testing

WARNING

Using a defective instrument can cause damage to the endoscope and result in injury to the patient or user.

NOTE

If the packaging is damaged before use, replace the instrument with a new one.

NOTE

Do not use damaged packaging or products that have expired.

NOTE

If you notice any irregularities, replace the instrument with a new one.

3.1.1 Compatibility and accessories

- **ZB0100:** OLYMPUS/ FUJIFILM endoscopes (color blue)
- **ZB0120:** PENTAX endoscopes (color red)

3.2. Application

1. Unpack the biopsy valve.
2. Securely attach the biopsy valve with the cap closed to the biopsy/suction channel opening of the endoscope.

3.3. Technical Data

The list of all variants as well as the materials and individual parts used can be found in the technical product data sheet and can be obtained on request from ENDO-FLEX GmbH the and/or on the website of the . ENDO-FLEX GmbH

4. Storage & Transportation

This product may only be transported and stored in the packaging provided for it. In addition, no other special transport conditions are required.

Store the product at room temperature and protect it from direct sunlight and moisture. Do not store instruments near chemicals, disinfectants, or radioactive radiation.

SAFETY INSTRUCTIONS

Do not place objects on the storage packaging and the sterile barrier system!

5. Disposal

After use, this product may present a biological hazard. Disposal must be carried out in accordance with the regulations for clinically contaminated materials and national and regional laws and regulations.

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