

# ReBite<sup>®</sup>

Biopsy Forceps - Reusable



# ReBite®

## Biopsy Forceps - Reusable

### Features:

- Detachable Cups
- Teflon® Coated Inner Pull Wire
- Ground Distal Spiral (10 cm)
- Autoclavable

### Benefits:

- Optimal Cleaning and Enhanced Instrument Life
- Smooth Operation of Cups
- Easy Passage through Working Channel
- Easy and Effective Sterilization

### Order Information:

Catalog No.	Description	Outer Diameter	Length	Minimum Working Channel Diameter
<b>Biopsy Forceps (with handle)</b>				
<b>Pulmonary</b>				
400220	Mini Biopsy Forceps, oval cups	1.0 mm	110 cm	1.2 mm
401240	Biopsy Forceps, oval cups	1.8 mm	120 cm	2.0 mm
402240	Biopsy Forceps, oval cups with needle	1.8 mm	120 cm	2.0 mm
405240	Biopsy Forceps, toothed cups	1.8 mm	120 cm	2.0 mm
406240	Biopsy Forceps, oval cups	2.3 mm	120 cm	2.5 mm
407240	Biopsy Forceps, oval cups with needle	2.3 mm	120 cm	2.5 mm
408240	Biopsy Forceps, toothed cups	2.3 mm	120 cm	2.5 mm
<b>Gastro &amp; Colon</b>				
400360	Mini Biopsy Forceps, oval cups	1.0 mm	180 cm	1.2 mm
400460	Mini Biopsy Forceps, oval cups	1.0 mm	230 cm	1.2 mm
401360	Biopsy Forceps, oval cups	1.8 mm	180 cm	2.0 mm
401460	Biopsy Forceps, oval cups	1.8 mm	230 cm	2.0 mm
402360	Biopsy Forceps, oval cups with needle	1.8 mm	180 cm	2.0 mm
402460	Biopsy Forceps, oval cups with needle	1.8 mm	230 cm	2.0 mm
405360	Biopsy Forceps, toothed cups	1.8 mm	180 cm	2.0 mm
405460	Biopsy Forceps, toothed cups	1.8 mm	230 cm	2.0 mm
406360	Biopsy Forceps, oval cups	2.3 mm	180 cm	2.5 mm
406460	Biopsy Forceps, oval cups	2.3 mm	230 cm	2.5 mm
407360	Biopsy Forceps, oval cups with needle	2.3 mm	180 cm	2.5 mm
407460	Biopsy Forceps, oval cups with needle	2.3 mm	230 cm	2.5 mm
408360	Biopsy Forceps, toothed cups	2.3 mm	180 cm	2.5 mm
408460	Biopsy Forceps, toothed cups	2.3 mm	230 cm	2.5 mm
408360WN	Biopsy Forceps, toothed cups with needle	2.3 mm	180 cm	2.5 mm
408460WN	Biopsy Forceps, toothed cups with needle	2.3 mm	230 cm	2.5 mm
<b>Jumbo – Forceps</b>				
411360	Jumbo - Biopsy Forceps, oval cups	3.4 mm	180 cm	3.6 mm
411460	Jumbo - Biopsy Forceps, oval cups	3.4 mm	230 cm	3.6 mm
412360	Jumbo - Biopsy Forceps, oval cups with needle	3.4 mm	180 cm	3.6 mm
412460	Jumbo - Biopsy Forceps, oval cups with needle	3.4 mm	230 cm	3.6 mm

# ReBite®

## Hot Biopsy Forceps - Reusable

### Features:

- Teflon® Coated Inner Pull Wire -----> Smooth Operation of Cups
- Teflon® Isolated Spiral -----> Easy Passage through Working Channel
- Handle and Forceps are Two Separate Components -> Permits Individual Component Replacement
- Completely Autoclavable -----> Easy and Effective Sterilization

### Benefits:


### Order Information:

Catalog No.	Description	Outer Diameter	Length	Minimum Working Channel Diameter
<b>Hot Biopsy Forceps (handle excluded)</b>				
<b>Gastro &amp; Colon</b>				
<b>406360H</b>	Biopsy Forceps, oval cups	2.5 mm	180 cm	2.7 mm
<b>406460H</b>	Biopsy Forceps, oval cups	2.5 mm	230 cm	2.7 mm
<b>407360H</b>	Biopsy Forceps, oval cups with needle	2.5 mm	180 cm	2.7 mm
<b>407460H</b>	Biopsy Forceps, oval cups with needle	2.5 mm	230 cm	2.7 mm

### Related Accessories required for Hot Biopsy Forceps:

Hot Biopsy Forceps are delivered without handle. When ordering a HF version, please add:

**Handle for all Hot Biopsy Forceps Catalog No. 100500**



Medi-Globe “Hot” handles are compatible with Olympus and Microvative electro-surgical generators. To ensure compatibility with other types of electro-surgical generators, please contact us or your local distributor. We recommend the use of compatible Medi-Globe active cords. You will find the corresponding Medi-Globe active cords under “Accessories”.

Manufactured by:

 **Medi-Globe**

Medi-Globe GmbH

Medi-Globe-Strasse 1-5

83101 Achenmühle, Germany

Phone: ++49 8032 973 379

Fax: ++49 8032 973 339

e-mail: [sales@medi-globe.de](mailto:sales@medi-globe.de)

[www.medi-globe.de](http://www.medi-globe.de)

Distributed by: