


Fibrin Needles

Single Use



Fibrin Needles

Single Use

Features:

- Double Lumen Coaxiale „Needle in Needle“-System
- Removable Safety Clip
- 15° Beveled Needle
- Easy Grip Handle
- Single Use

Benefits:

- Precise Confluence of Fibrin/Thrombin at Intended Position
- Secure Needle within the Catheter
- Effective Needle Puncture
- Easy Handling
- Safety and Convenience

Order Information:



Catalog No.	Needle Diameter	Needle Length	Catheter Diameter	Total Length
GFN-01-14-180	1.2 mm	4.0 mm	2.5 mm	180 cm
GFN-01-04-230	1.2 mm	4.0 mm	2.5 mm	230 cm
GFN-01-06-180	1.2 mm	6.0 mm	2.5 mm	180 cm
GFN-01-06-230	1.2 mm	6.0 mm	2.5 mm	230 cm
GFN-01-08-180	1.2 mm	8.0 mm	2.5 mm	120 cm
GFN-01-08-230	1.2 mm	8.0 mm	2.5 mm	180 cm

Manufactured by:

Medi-Globe
 Medi-Globe GmbH
 Medi-Globe-Strasse 1-5
 83101 Achenmühle, Germany
 Phone: ++49 8032 973 379
 Fax: ++49 8032 973 339
 e-mail: sales@medi-globe.de
 www.medi-globe.de

Distributed by: