

# Nasal Biliary Drainages

Single Use



# Nasal Biliary Drainages

Single Use


## Features:

- Various Types and Sizes Available
- Tapered and Rounded Ends
- Well-Balanced Wall Thickness
- Radiopaque Material
- Single Use

## Benefits:

- Fits all Anatomic Cases
- Atraumatic Implantation, Easy Placement
- Maximum Inner Lumen with Guaranteed Form Stability and Kink Resistance
- Exact Positioning
- Safety and Convenience

## Nasal Biliary Drainage Sets:



## Order Information:

### Nasal Biliary Drainage Sets

| Catalog No.           | Description           | Drainage Catheter |        | Nasal Tube | Drainage Connecting Tube Length | Recommended Guide Wire |        |
|-----------------------|-----------------------|-------------------|--------|------------|---------------------------------|------------------------|--------|
|                       |                       | French size       | Length | Length     |                                 | Diameter               | Length |
| <b>GND-01-05-290S</b> | Pigtail               | 5.0               | 290 cm | 50 cm      | 135 cm                          | .035" (0.89 mm)        | 400 cm |
| <b>GND-01-07-290S</b> | Pigtail               | 7.0               | 290 cm | 50 cm      | 135 cm                          | .035" (0.89 mm)        | 400 cm |
| <b>GND-01-10-290S</b> | Pigtail               | 10.0              | 290 cm | 50 cm      | 135 cm                          | .052" (1.32 mm)        | 400 cm |
| <b>GND-02-05-290S</b> | Standard              | 5.0               | 290 cm | 50 cm      | 135 cm                          | .035" (0.89 mm)        | 400 cm |
| <b>GND-02-07-290S</b> | Standard              | 7.0               | 290 cm | 50 cm      | 135 cm                          | .035" (0.89 mm)        | 400 cm |
| <b>GND-02-10-290S</b> | Standard              | 10.0              | 290 cm | 50 cm      | 135 cm                          | .052" (1.32 mm)        | 400 cm |
| <b>GND-03-05-290S</b> | Standard with Pigtail | 5.0               | 290 cm | 50 cm      | 135 cm                          | .035" (0.89 mm)        | 400 cm |
| <b>GND-03-07-290S</b> | Standard with Pigtail | 7.0               | 290 cm | 50 cm      | 135 cm                          | .035" (0.89 mm)        | 400 cm |
| <b>GND-03-10-290S</b> | Standard with Pigtail | 10.0              | 290 cm | 50 cm      | 135 cm                          | .052" (1.32 mm)        | 400 cm |

All Nasal Biliary Drainage Sets include a Drainage Catheter, Nasal Tube, Drainage Connecting Tube, Injection Connector and a Teflon® coated Guide Wire.

All Nasal Biliary Drainage Catheters are also available individually (catalog no. without „S“).



#### Nasal Tube

| Catalog No.          | French size | Length |
|----------------------|-------------|--------|
| <b>GND-41-05-050</b> | 5.0         | 50 cm  |
| <b>GND-41-07-050</b> | 7.0         | 50 cm  |
| <b>GND-41-10-050</b> | 10.0        | 50 cm  |




#### Injection Connector

| Catalog No.          | French size |
|----------------------|-------------|
| <b>GND-43-05-000</b> | 5.0         |
| <b>GND-43-07-000</b> | 7.0         |
| <b>GND-43-10-000</b> | 10.0        |



#### Guide Wire with Teflon®-coating, fixed core

| Catalog No.          | Diameter        | Length |
|----------------------|-----------------|--------|
| <b>GGW-01-35-400</b> | .035" (0.89 mm) | 400 cm |
| <b>GGW-01-52-400</b> | .052" (1.32 mm) | 400 cm |



#### Drainage Connecting Tube

| Catalog No.          | Length |
|----------------------|--------|
| <b>GND-42-03-135</b> | 135 cm |

Manufactured by:

 **Medi-Globe**

Medi-Globe GmbH

Medi-Globe-Strasse 1-5

83101 Achenmühle, Germany

Phone: ++49 8032 973 379

Fax: ++49 8032 973 339

e-mail: [sales@medi-globe.de](mailto:sales@medi-globe.de)

[www.medi-globe.de](http://www.medi-globe.de)

Distributed by: