




E.R.C.P. Catheters

Reusable



E.R.C.P. Catheters

Reusable


Features:

- Teflon® Tip with Memory-Effect
- Color Markings
- Metal Tips / Metal rings
- Removable Stylet (1.2 m)
- Autoclavable

Benefits:

- > Allows Pre-Shaping of the Tip
- > Macroscopically Controlled Cannulation Depth
- > Maximum X-Ray Visualization
- > Kink Free Cannulation
- > Easy and Effective Sterilization

Order Information:



Catalog No.	Description	Outer Diameter	Accepted Guide Wire
522420	Tapered Tip 10 mm	1.8 mm	.021" (0.53 mm)
522421	Tapered Tip 5 mm	1.8 mm	.021" (0.53 mm)
524420	Cone Metal Tip	1.8 mm	.021" (0.53 mm)
520420	Standard	1.8 mm	.035" (0.89 mm)
523420	Tapered Tip 10 mm	2.2 mm	.035" (0.89 mm)
523421	Tapered Tip 5 mm	2.2 mm	.035" (0.89 mm)
532420	"Bullet-Type", Metal Tip	1.8 mm	.035" (0.89 mm)
533420	"Ball Type", Metal Tip	1.8 mm	.035" (0.89 mm)
534420	Minor E.R.C.P. Catheter, Metal Cannula Length 1 mm	1.8 mm	N/A

Manufactured by:



Medi-Globe GmbH
 Medi-Globe-Strasse 1-5
 83101 Achenmühle, Germany
 Phone: ++49 8032 973 379
 Fax: ++49 8032 973 339
 e-mail: sales@medi-globe.de
 www.medi-globe.de

Distributed by: