

Premier Cut & S.O.S. Sphincterotomes

Single Use



**PREMIER CUT
SPHINCTEROTOME**
with SEMI-INSULATED CUTTING WIRE

S.O.S. SPHINCTEROTOME
with PRE-LOADED
GUIDE WIRE



Premier Cut & S.O.S. Sphincterotomes

Single Use

Features:

- 2 & 3 Lumen Versions
- Semi-Insulated Wire
- "Pre-curved" Versions
- Rotatable Cutting Wire
- Color Marking with Black Rings
- Ergonomically Shaped Handle

Benefits:

- > Contrast Medium Injection while Guide Wire is in Place
- > Safe and Controlled Papillary Roof Incision
- > Easier Positioning to the Papilla
- > Precise Positioning
- > Macroscopically Controlled Cannulation Depth
- > User-Friendly Handling

Order Information:

Premier Cut

3-Lumen Guide Wire Sphincterotome with semi-insulated cutting wire, rotatable

Catalog No.	Description	Outer Diameter	Cutting Wire Length	Guide Wire
GSP-27-20-030	Tapered Tip, Tip Length 5 mm	2.55 / 2.05 mm	30 mm	.035" (0.89 mm)

Premier Cut

2-Lumen Guide Wire Sphincterotome with semi-insulated cutting wire, rotatable

GSP-25-20-020	Tapered Tip, Tip Length 5 mm	2.55 / 2.05 mm	20 mm	.035" (0.89 mm)
GSP-25-20-030	Tapered Tip, Tip Length 5 mm	2.55 / 2.05 mm	30 mm	.035" (0.89 mm)
GSP-26-20-020	Tapered Tip, Tip Length 15 mm	2.55 / 2.05 mm	20 mm	.035" (0.89 mm)
GSP-26-20-030	Tapered Tip, Tip Length 15 mm	2.55 / 2.05 mm	30 mm	.035" (0.89 mm)

S.O.S. "Speed Option Set" Sphincterotome

2-Lumen Guide Wire Sphincterotome with pre-loaded hydrophilic Teflon® coated .035" (0.89 mm) nitinol guide wire, Length 450 cm (Catalog no. GGW-03-35-450)

Catalog No.	Description	Outer Diameter	Cutting Wire Length	Guide Wire
GSP-21-20-031	Tapered Tip, Tip Length 5 mm	2.55 / 2.05 mm	30 mm	Pre-loaded
GSP-25-20-031	Semi-insulated cutting wire, Tapered Tip, Tip Length 5 mm	2.55 / 2.05 mm	30 mm	Pre-loaded

Manufactured by:

Medi-Globe
 Medi-Globe GmbH
 Medi-Globe-Strasse 1-5
 83101 Achenmühle, Germany
 Phone: ++49 8032 973 379
 Fax: ++49 8032 973 339
 e-mail: sales@medi-globe.de
 www.medi-globe.de

Distributed by: