

# RotaCut<sup>®</sup> - M

Monofilament Sphincterotomes - Single Use

## MONOFILAMENT CUTTING WIRE

- > Precise & accurate cut
- > Smooth-running sphincterotomy



3-Lumen Version

 **Medi-Globe<sup>®</sup>**  
The Spirit of Care.™



# RotaCut<sup>®</sup> - M

Monofilament Sphincterotomes - Single Use

## Features:

## Benefits

- Monofilament Cutting Wire -----> Precise & Smooth Incision of the Papillary Roof
- 3-Lumen Version -----> Contrast Medium Injection while Guide Wire is in Place
- Rotatable Cutting Wire -----> Precise Positioning
- Color Markings with Black Rings -----> Macroscopically Controlled Cannulation Depth
- "Pre-curved" Versions -----> Easier Positioning to the Papilla
- Ergonomically Shaped Handle -----> User-Friendly Handling

## Order Information:

### 3-Lumen Guide Wire Sphincterotomes with monofilament cutting wire



Catalog No.	Description	Outer Diameter	Cutting Wire Length	Guide Wire
<b>GSP-32-20-020</b>	Tapered Tip, Tip Length 5 mm	2.55 / 2.05 mm	20 mm	.035" (0.89 mm)
<b>GSP-32-20-025</b>	Tapered Tip, Tip Length 5 mm	2.55 / 2.05 mm	25 mm	.035" (0.89 mm)
<b>GSP-32-20-030</b>	Tapered Tip, Tip Length 5 mm	2.55 / 2.05 mm	30 mm	.035" (0.89 mm)

Manufactured by:

 **Medi-Globe**

Medi-Globe GmbH  
Medi-Globe-Strasse 1-5  
83101 Achenmühle, Germany  
Phone: ++49 8032 973 379  
Fax: ++49 8032 973 339  
e-mail: sales@medi-globe.de  
www.medi-globe.de

Distributed by: