

Stone Extraction Baskets

Reusable



Stone Extraction Baskets

Reusable

Features:

- Adjustable Sealing Cap
- Ergonomic Instrument Holder
- All Components are Autoclavable

Benefits:


- Contrast Medium Compatibility
- Simple and Safe Manipulation of the Basket
- Easy and Effective Sterilization

Order Information:

Stone Extraction Basket Set – consists of Basket, Teflon® Tube, Instrument Holder and Sealing Cap

| Catalog No. | Description | Outer Diameter | Basket Diameter |
|--------------------------------|---------------|----------------|-----------------|
| Basket with Tube 180 cm | | | |
| 606420 | 4-Wire Basket | 1.9 mm | 20 mm |
| 606425 | 4-Wire Basket | 1.9 mm | 25 mm |
| 606430 | 4-Wire Basket | 1.9 mm | 30 mm |
| 608620 | 6-Wire Basket | 1.9 mm | 20 mm |
| 608625 | 6-Wire Basket | 1.9 mm | 25 mm |
| 608630 | 6-Wire Basket | 1.9 mm | 30 mm |
| 605420 | 4-Wire Basket | 2.5 mm | 20 mm |
| 605425 | 4-Wire Basket | 2.5 mm | 25 mm |
| 605430 | 4-Wire Basket | 2.5 mm | 30 mm |
| 607620 | 6-Wire Basket | 2.5 mm | 20 mm |
| 607625 | 6-Wire Basket | 2.5 mm | 25 mm |
| 607630 | 6-Wire Basket | 2.5 mm | 30 mm |
| Basket with Tube 280 cm | | | |
| 615420 | 4-Wire Basket | 2.5 mm | 20 mm |
| 615425 | 4-Wire Basket | 2.5 mm | 25 mm |
| 615430 | 4-Wire Basket | 2.5 mm | 30 mm |
| 615620 | 6-Wire Basket | 2.5 mm | 20 mm |
| 615625 | 6-Wire Basket | 2.5 mm | 25 mm |

As an alternative to the instrument holder the handle order no. 100600 can be used.



Manufactured by:

Medi-Globe

Medi-Globe GmbH

Medi-Globe-Strasse 1-5

83101 Achenmühle, Germany

Phone: ++49 8032 973 379

Fax: ++49 8032 973 339

e-mail: sales@medi-globe.de

www.medi-globe.de

Distributed by: