

# RotaCut<sup>®</sup>

Rotatable Sphincterotomes - Reusable



# RotaCut®

Rotatable Sphincterotomes - Reusable

## Features:

- Rotatable Catheter Tip
- Double Lumen Versions
- Color Markings with Black Rings
- Detachable Catheter and Handle

## Benefits:

- Precise Positioning
- Guide Wire or Contrast Medium Compatibility
- Macroscopically Controlled Cannulation Depth and Cutting Wire Position
- Economical Component Replacement

## Order Information:

### 2-Lumen Guide Wire Sphincterotomes, rotatable:



Catalog No.	Description	Outer Diameter	Cutting Wire Length	Guide Wire
512420	Tapered Tip, Tip Length 5 mm	2.55/2.05 mm	20 mm	.035"
512430	Tapered Tip, Tip Length 5 mm	2.55/2.05 mm	30 mm	.035"
519420	B II Variable Cutting Wire	2.55/2.05 mm	0-20 mm	.035"

### Standard Sphincterotomes, 1-Lumen, rotatable:



Catalog No.	Description	Outer Diameter	Cutting Wire Length	Tip Length
502420	Standard	1.7 mm	20 mm	5 mm
502430	Standard	1.7 mm	30 mm	5 mm
503420	Tapered Tip	1.7 mm	20 mm	5 mm
503430	Tapered Tip	1.7 mm	30 mm	5 mm

### Special Version, non-rotatable

Catalog No.	Description	Outer Diameter	Cutting Wire Length
516415	Needle Knife Sphincterotome, variable Needle Length	1.7 mm	0-15 mm

### Required Handle

Catalog No.	Description
100520	Handle, for Sphincterotomes – Reusable

Medi-Globe handles are compatible with Olympus electro-surgical generators. To ensure compatibility with other types of electro-surgical generators, please contact us or your local distributor.

We recommend the use of compatible Medi-Globe active cords. You will find the corresponding Medi-Globe active cords in the leaflet "Accessories".

Manufactured by:

**Medi-Globe®**

Medi-Globe GmbH  
 Medi-Globe-Strasse 1-5  
 83101 Achenmühle, Germany  
 Phone: +49 8032 973 379  
 Fax: +49 8032 973 339  
 e-mail: sales@medi-globe.de  
 www.medi-globe.de

Distributed by: