

Accessories

Spray Catheters & Active Cords

Spray Catheters

Application of Colorizing Liquids for
Chromoendoscopy & Narcotics.

Active Cords

Connecting "Hot"-Biopsy Forceps,
Polypectomy Snares and Sphincterotomes to
HF-Generators of Various Manufacturers.



Accessories

Spray Catheters & Active Cords

Spray Catheters

Coaxial Spray Catheter for Chromoendoscopy and the Application of Narcotics.

Spray Catheter – Single Use

Catalog No.	Spray Head Diameter	Catheter Diameter	Working Length
GSC-01-25-180	2.5 mm	1.8 mm	180 cm
GSC-01-25-230	2.5 mm	1.8 mm	230 cm

Spray Catheter – Reusable

Catalog No.	Spray Head Diameter	Catheter Diameter	Working Length
540360	2.3 mm	1.8 mm	180 cm
540460	2.3 mm	1.8 mm	230 cm

Active Cords

Active Cords for Connecting Medi-Globe Instruments to HF-Generators of Various Manufacturers.

Active Cord

Catalog No.	Description	Length
GAC-51-04-300	Active Cord, Erbe 4 mm	300 cm
GAC-51-05-300	Active Cord, Erbe 5 mm	300 cm
GAC-52-04-300	Active Cord, Martin / Berchthold 4 mm	300 cm
GAC-53-08-300	Active Cord, Valleylab 8 mm	300 cm
GAC-54-06-300	Active Cord, Olympus 6 mm	300 cm

Distributed by: