

FlexiBite

Biopsy Forceps – Single Use

FlexiBite

- Sharp Standard & Serrated Cups
Precise & Clean Tissue Acquisition
- Spiked Versions
Controlled & Accurate Sampling from the Targeted Biopsy Site
- Versions with Coated Metal Spiral with Positioning Bands
Friction Minimized Passage through Working Channel
- Convenient Evaluation of Forceps Position for a Controlled & Safe Catheter Advancement and Removal*
- Hot-Biopsy Forceps
Safe Electro-coagulation and Excision of Small Polyps



Cup Shapes



Oval Cups



Oval Cups with Spike



Serrated Oval Cups



Serrated Oval Cups with Spike

Biopsy Forceps Standard



Oval Cups



Oval Cups with Spike

FlexiBite – Biopsy Forceps with Metal Spiral, sterile packed

Catalog No.	Cup Shapes	Outer Diameter	Working Length	Minimum Working Channel Diameter
GBF-01-18-120	Oval Cups	1.8 mm	120 cm	2.0 mm
GBF-01-23-180	Oval Cups	2.3 mm	180 cm	2.5 mm
GBF-01-23-230	Oval Cups	2.3 mm	230 cm	2.5 mm
GBF-02-23-180	Oval Cups with Spike	2.3 mm	180 cm	2.5 mm
GBF-02-23-230	Oval Cups with Spike	2.3 mm	230 cm	2.5 mm



Serrated Oval Cups



Serrated Oval Cups with Spike

FlexiBite – Biopsy Forceps with Metal Spiral, sterile packed

Catalog No.	Cup Shapes	Outer Diameter	Working Length	Minimum Working Channel Diameter
GBF-03-23-230	Serrated Oval Cups	2.3 mm	230 cm	2.5 mm
GBF-04-23-230	Serrated Oval Cups with Spike	2.3 mm	230 cm	2.5 mm

Biopsy Forceps with Coated Metal Spiral



FlexiBite – Biopsy Forceps with Coated Metal Spiral with Positioning Bands, sterile packed

Catalog No.	Cup Shapes	Outer Diameter	Working Length	Minimum Working Channel Diameter
GBF-21-23-180	Oval Cups	2.3 mm	180 cm	2.5 mm
GBF-21-23-230	Oval Cups	2.3 mm	230 cm	2.5 mm
GBF-22-23-180	Oval Cups with Spike	2.3 mm	180 cm	2.5 mm
GBF-22-23-230	Oval Cups with Spike	2.3 mm	230 cm	2.5 mm



FlexiBite – Biopsy Forceps with Coated Metal Spiral with Positioning Bands, sterile packed

Catalog No.	Cup Shapes	Outer Diameter	Working Length	Minimum Working Channel Diameter
GBF-23-23-230	Serrated Oval Cups	2.3 mm	230 cm	2.5 mm
GBF-24-23-230	Serrated Oval Cups with Spike	2.3 mm	230 cm	2.5 mm



Positioning Bands on the Coating

“Hot”-Biopsy Forceps




Monopolar “Hot”-Biopsy Forceps with Coated Metal Spiral, sterile packed

Catalog No.	Cup Shapes	Outer Diameter	Working Length	Minimum Working Channel Diameter
GBF-11-24-180	Oval Cups	2.4 mm	180 cm	2.8 mm
GBF-11-24-230	Oval Cups	2.4 mm	230 cm	2.8 mm

Medi-Globe “Hot”-Biopsy Forceps are directly compatible with Olympus electro-surgical generators. To verify compatibility with other electro-surgical generators please contact your Medi-Globe sales representative. Please also consult the Medi-Globe “Accessories” leaflet to see a complete list of active cords manufactured to fit most generators.

FlexiBite

Biopsy Forceps – Single Use



Distributed by:

Ver. 4.0 (2010-06)

Manufactured by: