

Endoscopic Ultrasound FNA Systems

For Pentax Fibre Ultrasound Endoscopes FG-32UA, FG-34UX, FG-36UX & FG-38UX

Endoscopic Ultrasound “Hancke-Vilmann” FNA Systems

- Reusable Components (Handle & Metal Spiral)
Combined with Single Use Needle / Stylet Unit
Cost-effective
Individual Component Replacement
- Echogenic Needle Tip Preparation
Highly Visible Ultrasound Imaging
- Graduated Needle Adjuster
Variable Needle Projection (0 - 12 cm)



Endoscopic Ultrasound FNA Systems

Endoscopic Ultrasound FNA Systems

These EUS-FNA systems, also well known as type “Hancke-Vilmann” needles, are particularly designed and manufactured for a specific range of Pentax Fibre Ultrasound Endoscopes. A complete EUS-FNA “Hancke-Vilmann” needle set consist of an aluminum handle and a metal spiral, which are reusable and autoclavable and a disposable aspiration needle / stylet unit, which has to discard after the treatment. All components are supplied sterile separately and have to assemble prior of use.

22 Gauge Needles

EUS-FNA “Hancke-Vilmann” Needle Set for Pentax FG-32UA, sterile packed

Catalog No.	Description	Metal Spiral Outer Diameter	Minimum Accessory Channel Diameter
222280	Set includes: 1 x Handle (100202), Reusable 1 x Metal Spiral (222284), Reusable 1 x Aspiration Needle / Stylet Unit (222282), Single Use	1.8 mm	2.0 mm

Replacement Components for EUS-FNA “Hancke-Vilmann” Needle Set for Pentax FG-32UA

Catalog No.	Description	Metal Spiral Outer Diameter	Minimum Accessory Channel Diameter
222282	Aspiration Needle / Stylet Unit, Single Use, sterile packed	-	-
222284	Metal Spiral, Reusable, sterile packed	1.8 mm	2.0 mm

EUS-FNA “Hancke-Vilmann” Needle Set for Pentax FG-34UX and FG-36UX, sterile packed

Catalog No.	Description	Metal Spiral Outer Diameter	Minimum Accessory Channel Diameter
222480	Set includes: 1 x Handle (100202), Reusable 1 x Metal Spiral (222484), Reusable 1 x Aspiration Needle / Stylet Unit (222482), Single Use	1.8 mm	2.0 mm

Replacement Components for EUS-FNA “Hancke-Vilmann” Needle Set for Pentax FG-34UX and FG-36UX

Catalog No.	Description	Metal Spiral Outer Diameter	Minimum Accessory Channel Diameter
222482	Aspiration Needle / Stylet Unit, Single Use, sterile packed	-	-
222484	Metal Spiral, Reusable, sterile packed	1.8 mm	2.0 mm

19 Gauge Needle

EUS-FNA “Hancke-Vilmann” Needle Set for Pentax FG-38UX, sterile packed

Catalog No.	Description	Metal Spiral Outer Diameter	Minimum Accessory Channel Diameter
223280	Set includes: 1 x Handle (100202), Reusable 1 x Metal Spiral (223284), Reusable 1 x Aspiration Needle / Stylet Unit (223282), Single Use	2.1 mm	2.3 mm

Replacement Components for EUS-FNA “Hancke-Vilmann” Needle Set for Pentax FG-38UX


Catalog No.	Description	Metal Spiral Outer Diameter	Minimum Accessory Channel Diameter
223282	Aspiration Needle / Stylet Unit, Single Use, sterile packed	-	-
223484	Metal Spiral, Reusable, sterile packed	2.1 mm	2.3 mm

Handle for EUS-FNA “Hancke-Vilmann” Needle Systems


Catalog No.	Description
100202	Aluminium Handle with Graduated Needle Adjuster, Reusable, sterile packed

Accessories for 19 Gauge Endoscopic Ultrasound FNA Systems

Nitinol Guide Wire Compatible to 19 Gauge EUS-FNA Systems


Catalog No.	Description	Diameter	Length	Nitinol Guide Wire
GGW-08-30-400	Nitinol Guide Wire, uncoated with 15 cm Radiopaque Flexible Spring Tip	.030"	400 cm	

Pseudo-Cyst Drainage Stents, Double Pigtail with Position Markers

Catalog No.	Description	Size	Length	Pseudo-Cyst Drainage Stent
GBS-11-08-002	Double Pigtail Stent with Position Markers	8.5 FR	2 cm	
GBS-11-08-005		8.5 FR	5 cm	
GBS-11-08-007		8.5 FR	7 cm	
GBS-11-08-009		8.5 FR	9 cm	
GBS-11-10-002		10.0 FR	2 cm	
GBS-11-10-005		10.0 FR	5 cm	
GBS-11-10-007		10.0 FR	7 cm	
GBS-11-10-009		10.0 FR	9 cm	

Endoscopic Ultrasound FNA Systems

For Pentax Fibre Ultrasound Endoscopes FG-32UA, FG-34UX, FG-36UX & FG-38UX



Distributed by:

Ver. 3.0 (2010-10)

Manufactured by: