

Aspiration Needles

Single Use



Safety Clip and Locking Mechanism

Secure needle advancement and protection

15° Beveled Needle with Side Hole

Atraumatic piercing capability for maximum yield of aspirates

Aspiration Needles – Single Use



Aspiration Needle for Bronchoscopy

With metal spiral; including 10ml vacuum syringe with stopcock

Catalogue No.	Needle Diameter (mm/Gauge)	Needle Length (mm)	Metal Spiral Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GAN-01-19-100	0.8/21	12	1.9	100	2.1



L-GAN-SU-00-EN / LL-0009 Vers. 6.0 (12/2020)



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
Email: sales@medi-globe.de • www.medi-globe.de