

Multiband-Ligator

For Esophageal Varices – Single Use
Non Sterile



User-friendly Mounting

Device fixation onto endoscopes without adapter

Rapid & Easy Operation of the “Shooting” Mechanism

Facilitates band release in emergency situations

Colored Alert Band

Indicates second to last band

Multiband-Ligator – Single Use



Esophageal Multiband Ligator

Pre-loaded barrel with 5 bands (contains latex); improved mounting to the endoscope; additional Luer-Lock irrigation port; maximum endoscope length 150 cm

Catalogue No.	Suitable for Endoscope Outer Diameter (mm)	Number of Bands	Band Inner Diameter in Released Position (mm)	Length of the Introduction Catheter (cm)	Min. Accessory Channel (mm)
GBL-01-05-001	8.5–11.5	5	1.9	150	2.8

L-GBL-SU-00-EN / LL-0238 Vers. 2.0 (12/2020)



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
Email: sales@medi-globe.de • www.medi-globe.de

