

Pancreatic Stents – Single Components – Single Use

Pancreatic Stents and Stent-Sets


Single Use



Pancreatic Stent – Straight with Single Pigtail

Single pigtail without flap; multiple drainage holes; 1 position marker

Catalogue No.	Size(French)	Length between Pigtail (cm)
GBS-18-05-003	5.0	3
GBS-18-05-004	5.0	4
GBS-18-05-005	5.0	5
GBS-18-05-006	5.0	6
GBS-18-05-007	5.0	7



Pancreatic Stent – S-shaped

Single components

Catalogue No.	Size (French)	Length (cm)
GBS-09-05-005	5.0	5
GBS-09-05-007	5.0	7
GBS-09-05-009	5.0	9
GBS-09-05-012	5.0	12
GBS-09-05-015	5.0	15
GBS-09-07-005	7.0	5
GBS-09-07-007	7.0	7
GBS-09-07-009	7.0	9
GBS-09-07-012	7.0	12
GBS-09-07-015	7.0	15

L-GBS-SU-01-EN / LL-0019 Vers. 6.0 (01/2021)



Balanced Wall Diameters

Maximum inner lumen with form stability and kink resistance

Tapered and Round Endings

Atraumatic implantation

Versions with Multiple Holes


For drainage of obstructed pancreatic ducts



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
Email: sales@medi-globe.de • www.medi-globe.de

Pancreatic Stent Sets – Single Use




Pancreatic Stent-Set – Straight

Straight with radial holes; 1 distal and 1 duodenal flap; compatible with guide wire .035"

Catalogue No.	Stent		Pusher	
	Size (French)	Length between Flaps (cm)	Size (French)	Length (cm)
GBS-10-05-003S	5.0	3	5.0	220
GBS-10-05-005S	5.0	5	5.0	220
GBS-10-05-007S	5.0	7	5.0	220
GBS-10-05-009S	5.0	9	5.0	220
GBS-10-05-011S	5.0	11	5.0	220
GBS-10-07-003S	7.0	3	7.0	220
GBS-10-07-005S	7.0	5	7.0	220
GBS-10-07-007S	7.0	7	7.0	220
GBS-10-07-008S	7.0	8	7.0	220
GBS-10-07-009S	7.0	9	7.0	220
GBS-10-07-011S	7.0	11	7.0	220
GBS-10-07-016S	7.0	16	7.0	220
GBS-10-08-005S	8.5	5	8.5	220
GBS-10-08-007S	8.5	7	8.5	220
GBS-10-08-008S	8.5	8	8.5	220
GBS-10-08-009S	8.5	9	8.5	220
GBS-10-08-011S	8.5	11	8.5	220
GBS-10-08-016S	8.5	16	8.5	220


Pancreatic Stent Sets – Single Use



Pancreatic Stent-Set – Straight

Straight with radial holes; 1 duodenal flap; compatible with guide wire .035"

Catalogue No.	Stent		Pusher	
	Size (French)	Length between Flaps (cm)	Size (French)	Length (cm)
GBS-20-05-003S	5.0	3	5.0	220
GBS-20-05-005S	5.0	5	5.0	220
GBS-20-05-007S	5.0	7	5.0	220



Pancreatic Stent-Set – S-shaped

Compatible with guide wire .035"

Catalogue No.	Stent		Pusher	
	Size (French)	Length (cm)	Size (French)	Length (cm)
GBS-09-05-005S	5.0	5	5.0	220
GBS-09-05-007S	5.0	7	5.0	220
GBS-09-05-009S	5.0	9	5.0	220
GBS-09-05-012S	5.0	12	5.0	220
GBS-09-05-015S	5.0	15	5.0	220
GBS-09-07-005S	7.0	5	7.0	220
GBS-09-07-007S	7.0	7	7.0	220
GBS-09-07-009S	7.0	9	7.0	220
GBS-09-07-012S	7.0	12	7.0	220
GBS-09-07-015S	7.0	15	7.0	220

Pancreatic Stent Sets – Single Use



Pancreatic Stent-Set – Straight with Single Pigtail

Single pigtail with multiple flaps; multiple drainage holes; 1 position marker; compatible with guide wire .035"

Catalogue No.	Stent		Pusher	
	Size (French)	Length between Pigtail (cm)	Size (French)	Length (cm)
GBS-08-05-003S	5.0	3	5.0	220
GBS-08-05-004S	5.0	4	5.0	220
GBS-08-05-005S	5.0	5	5.0	220
GBS-08-05-006S	5.0	6	5.0	220
GBS-08-05-007S	5.0	7	5.0	220
GBS-08-05-008S	5.0	8	5.0	220
GBS-08-05-009S	5.0	9	5.0	220
GBS-08-05-010S	5.0	10	5.0	220
GBS-08-05-012S	5.0	12	5.0	220
GBS-08-07-002S	7.0	2	7.0	220
GBS-08-07-004S	7.0	4	7.0	220
GBS-08-07-006S	7.0	6	7.0	220
GBS-08-07-009S	7.0	9	7.0	220
GBS-08-07-012S	7.0	12	7.0	220
GBS-08-08-006S	8.5	6	8.5	220
GBS-08-08-009S	8.5	9	8.5	220
GBS-08-08-012S	8.5	12	8.5	220
GBS-08-10-006S	10.0	6	10.0	220
GBS-08-10-009S	10.0	9	10.0	220
GBS-08-10-012S	10.0	12	10.0	220
GBS-08-11-006S	11.5	6	11.5	220
GBS-08-11-009S	11.5	9	11.5	220
GBS-08-11-012S	11.5	12	11.5	220

Pancreatic Stent Sets and Single Components – Single Use



Pancreatic Stent-Set – Straight with Single Pigtail

Without flap; multiple drainage holes; 1 position marker; compatible with guide wire .035"

Catalogue No.	Stent		Pusher	
	Size (French)	Length between Pigtail (cm)	Size (French)	Length (cm)
GBS-18-05-003S	5.0	3	5.0	220
GBS-18-05-004S	5.0	4	5.0	220
GBS-18-05-005S	5.0	5	5.0	220
GBS-18-05-006S	5.0	6	5.0	220
GBS-18-05-007S	5.0	7	5.0	220




Pusher

Single components; 1 metal ring

Catalogue-No.	Size (French)	Color Code	Length (cm)
GBS-42-05-220	5.0	grey	220
GBS-42-07-160	7.0	yellow	160
GBS-42-07-220	7.0	yellow	220
GBS-42-08-160	8.5	orange	160
GBS-42-08-180	8.5	orange	180
GBS-42-08-220	8.5	orange	220
GBS-42-10-160	10.0	red	160
GBS-42-10-180	10.0	red	180
GBS-42-10-220	10.0	red	220
GBS-42-11-160	11.5	blue	160
GBS-42-11-180	11.5	blue	180
GBS-42-11-220	11.5	blue	220

Pancreatic Stents – Single Components – Single Use



Pancreatic Stent – Straight, 2 Flaps

With radial holes; 1 distal and 1 duodenal flap

Catalogue No.	Size (French)	Length between Flaps (cm)
GBS-10-05-003	5.0	3
GBS-10-05-005	5.0	5
GBS-10-05-007	5.0	7
GBS-10-05-009	5.0	9
GBS-10-05-011	5.0	11
GBS-10-07-003	7.0	3
GBS-10-07-005	7.0	5
GBS-10-07-007	7.0	7
GBS-10-07-008	7.0	8
GBS-10-07-009	7.0	9
GBS-10-07-011	7.0	11
GBS-10-07-016	7.0	16
GBS-10-08-005	8.5	5
GBS-10-08-007	8.5	7
GBS-10-08-008	8.5	8
GBS-10-08-009	8.5	9
GBS-10-08-011	8.5	11
GBS-10-08-016	8.5	16

Pancreatic Stent – Straight, 1 Flap

With radial holes; 1 duodenal flap

Catalogue No.	Size (French)	Length (cm)
GBS-20-05-003	5.0	3
GBS-20-05-005	5.0	5
GBS-20-05-007	5.0	7

Pancreatic Stents – Single Components – Single Use



Pancreatic Stent – Straight with Single Pigtail

Single pigtail with multiple flaps; multiple drainage holes; 1 position marker

Catalogue No.	Size (French)	Length between Pigtail (cm)
GBS-08-05-003	5.0	3
GBS-08-05-004	5.0	4
GBS-08-05-005	5.0	5
GBS-08-05-006	5.0	6
GBS-08-05-007	5.0	7
GBS-08-05-008	5.0	8
GBS-08-05-009	5.0	9
GBS-08-05-010	5.0	10
GBS-08-05-012	5.0	12
GBS-08-07-002	7.0	2
GBS-08-07-004	7.0	4
GBS-08-07-006	7.0	6
GBS-08-07-009	7.0	9
GBS-08-07-012	7.0	12
GBS-08-08-006	8.5	6
GBS-08-08-009	8.5	9
GBS-08-08-012	8.5	12
GBS-08-10-006	10.0	6
GBS-08-10-009	10.0	9
GBS-08-10-012	10.0	12
GBS-08-11-006	11.5	6
GBS-08-11-009	11.5	9
GBS-08-11-012	11.5	12