



Medi-Globe®

A MEDI-GLOBE COMPANY

Crown Clip

Single Use



Synchronized Handle <> Jaw Rotation

For maneuverability and device orientation

Short Rigid and Compact Clip Head

Facilitates a safe clip fixation

Performance Jaw Tips in “Crown” Design

Secure vessel ligation and tissue closure

Repositionable Jaw

Allows clip adjustment prior to final release

Crown Clip – Single Use



Gastroscopy

Hemostasis clip system with coated metal spiral and with positioning bands*

Catalogue No.	Clip Opening (mm)	Working Length (cm)	Minimum Accessory Channel (mm)
GCS-12-26-165	12	165	2.8




Colonoscopy

Hemostasis clip system with coated metal spiral and with positioning bands*

Catalogue No.	Clip Opening (mm)	Working Length (cm)	Minimum Accessory Channel (mm)
GCS-12-26-230	12	230	2.8



* Medi-Globe Crown Clips are classified as  conditional. For more detailed information, please regard the information for use and the corresponding patient card or visit our website www.medi-globe.de



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
Email: sales@medi-globe.de • www.medi-globe.de