


# ExPander Multi-Stage Dilation Balloons

Single Use



L-GDB-SU-01-EN / LL-0213 Vers. 3.0 (12/2020)

## Multi-Stage Dilation Concept

Balloon for gradual adjustment

## Soft Balloon Material and Balloon Tip – Latex free

Smooth and atraumatic insertion

## Radiopaque Markers

Visibility and position control under X-ray vision

## Available with Pre-loaded Guide Wire

Time-saving balloon guidance and advancement



### Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle  
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392  
Email: sales@medi-globe.de • www.medi-globe.de

# Multi-Stage Dilation Balloons – Single Use



## Esophageal Dilation Balloon – Guide Wire Compatibility

Pre-loaded hydrophilic guide wire .035"; with radiopaque markers; 3 progressive stages; total length 180 cm

Catalogue No.	Balloon Diameter (mm)	Balloon Length (mm)	Catheter Outer Diameter (French)	Max. Pressure (psi/atm)	Min. Accessory Channel (mm)
GDB-04-08-080	6-7-8	80	2.3	147 / 10	2.8
GDB-04-10-080	8-9-10	80	2.3	132 / 9	2.8
GDB-04-12-080	10-11-12	80	2.3	118 / 8	2.8
GDB-04-15-080	12-13.5-15	80	2.3	118 / 8	2.8
GDB-04-18-080	15-16.5-18	80	2.3	103 / 7	2.8
GDB-04-20-080	18-19-20	80	2.3	88 / 6	3.2



## Esophageal/Colonic Balloon – Guide Wire Compatibility

Pre-loaded hydrophilic guide wire .035"; with radiopaque markers; 3 progressive stages; total length 230 cm

Catalogue No.	Balloon Diameter (mm)	Balloon Length (mm)	Catheter Outer Diameter (French)	Max. Pressure (psi/atm)	Min. Accessory Channel (mm)
GDB-04-08-055	6-7-8	55	2.3	147 / 10	2.8
GDB-04-10-055	8-9-10	55	2.3	132 / 9	2.8
GDB-04-12-055	10-11-12	55	2.3	118 / 8	2.8
GDB-04-15-055	12-13.5-15	55	2.3	118 / 8	2.8
GDB-04-18-055	15-16.5-18	55	2.3	103 / 7	2.8
GDB-04-20-055	18-19-20	55	2.3	88 / 6	3.2

# Multi-Stage Dilation Balloons – Single Use



## Biliary Dilation Balloon – Guide Wire Compatibility

Pre-loaded hydrophilic guide wire .035"; with radiopaque markers; 3 progressive stages; total length 200 cm

Catalogue No.	Balloon Diameter (mm)	Balloon Length (mm)	Catheter Outer Diameter (French)	Max. Pressure (psi/atm)	Min. Accessory Channel (mm)
GDB-04-08-030	6-7-8	30	2.3	147 / 10	2.8
GDB-04-10-030	8-9-10	30	2.3	132 / 9	2.8
GDB-04-12-030	10-11-12	30	2.3	118 / 8	2.8
GDB-04-15-030	12-13.5-15	30	2.3	118 / 8	2.8



## Inflation/Deflation System\*

Filling medium: liquid

Catalogue No.	Description
FXB-70-20	Filling volume 70 ml, incl. manometer 0-20 atm, automatic plunger locking system allows comfortable operation, rotating Luer-Lock connector, transparent syringe barrel for easy detection and removal of air bubbles

\* This product might not be available in all countries. For more information please ask your local representative.