

Fibrin Needles

Single Use



**Double Lumen Coaxial
„Needle in Needle”-System**

Precise confluence of fibrinogen/thrombin
at intended position

Red Safety Clip

To secure needle within the catheter

Easy Grip Handle

For easy handling

Fibrin Needles – Single Use



Gastroscopy

Double lumen; coaxial needle; mixing of the adhesive components in the tissue

| Catalogue No. | Needle Diameter (mm) | Needle Length (mm) | Catheter Outer Diameter (mm) | Working Length (cm) | Min. Accessory Channel (mm) |
|---------------|----------------------|--------------------|------------------------------|---------------------|-----------------------------|
| GFN-01-04-180 | 1.2 | 4 | 2.5 | 180 | 2.8 |
| GFN-01-06-180 | 1.2 | 6 | 2.5 | 180 | 2.8 |



Colonoscopy

Double lumen; coaxial needle; mixing of the adhesive components in the tissue

| Catalogue No. | Needle Diameter (mm) | Needle Length (mm) | Catheter Outer Diameter (mm) | Working Length (cm) | Min. Accessory Channel (mm) |
|---------------|----------------------|--------------------|------------------------------|---------------------|-----------------------------|
| GFN-01-04-230 | 1.2 | 4 | 2.5 | 230 | 2.8 |
| GFN-01-06-230 | 1.2 | 6 | 2.5 | 230 | 2.8 |

L-GFN-SUJ-00-EN / LL-0025 Vers. 6.0 (02/2021)



Medi-Globe GmbH

Medi-Globe-Straße 1 –5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 - 973-0 • Email: sales@medi-globe.de
www.medi-globe.de

