



Cyst Guide Wire - Nitinol

Uncoated; compatible with 19 Gauge EUS FNA-Systems for pseudo-cyst drainage;
15 cm radiopaque and flexible tip

Catalogue-No.	Description	Diameter (inch/mm)	Length (cm)	Flexible Part (cm)	Colour
GGW-08-30-400	Radiopaque Tip	.030/0.75	400	15	Black



Standard Guide Wire

Stainless steel wire; Teflon-coated

Catalogue-No.	Description	Diameter (inch/mm)	Length (cm)	Flexible Part (cm)	Colour
GGW-01-35-450	Fixed Core	.035/0.89	450	n.a.	Green

L-GGW-SU-01-EN / LL-0024 Vers. 12.0 (12/2020)

Guide Wires

Single Use



Various Guide Wire Versions Available

- EasyPass Nitinol Guidewires
- Bavarian Nitinol Guidewires
- EUS Cyst Guidewires



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
Email: sales@medi-globe.de • www.medi-globe.de

EasyPass Nitinol Guide Wires – Single Use



EasyPass Nitinol Guide Wire

2-coloured; Teflon-coated; 5 cm hydrophilic tip; super-flexible; radiopaque

Catalogue-No.	Description	Diameter (inch/mm)	Length (cm)	Flexible Part (cm)	Colour
GGW-10-35-260	Straight Tip	.035/0.89	260	5	Blue/White
GGW-10-35-260J	Bended "J-Tip"	.035/0.89	260	5	Blue/White
GGW-10-35-450	Straight Tip	.035/0.89	450	5	Blue/White
GGW-10-35-450J	Bended "J-Tip"	.035/0.89	450	5	Blue/White

Bavarian Nitinol Guide Wires – Single Use



Bavarian Nitinol Guide Wire

2-coloured; hydrophilic and Teflon-coated; 5 cm radiopaque tip

Catalogue-No.	Description	Diameter (inch/mm)	Length (cm)	Flexible Part (cm)	Colour
GGW-03-21-400	Soft Tip	.021/0.53	400	5	Blue/White
GGW-03-21-450	Soft Tip	.021/0.53	450	5	Blue/White
GGW-03-25-450	Soft Tip	.025/0.64	450	5	Blue/White
GGW-03-35-400	Soft Tip	.035/0.89	400	5	Blue/White
GGW-03-35-450	Soft Tip	.035/0.89	450	5	Blue/White
GGW-04-35-400	Standard Tip	.035/0.89	400	5	Blue/White
GGW-04-35-450	Standard Tip	.035/0.89	450	5	Blue/White

Bavarian Nitinol Guide Wire – compatible with ShortWire Systems

2-coloured; hydrophilic and Teflon-coated; 5 cm radiopaque tip

Catalogue-No.	Description	Diameter (inch/mm)	Length (cm)	Flexible Part (cm)	Colour
GGW-03-25-260	Soft Tip	.025/0.64	260	5	Blue/White
GGW-03-35-260	Soft Tip	.035/0.89	260	5	Blue/White