

# SonoTip® ProControl

## EUS Fine Needle Aspiration (FNA) System – Single Use

L-GUS-SU-01-EN / LL-0108 Vers. 9.0 (02/2021)



### Twist-Lock Technology

For precise tube and needle length adjustment

### Laser Treated Needle Surface

For improved needle visibility

### Distal Needle Stabilization

Precise tube guiding through large working channels  
Stable needle position on the target tissue



### Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle  
Tel.: +49 (0) 8032 - 973 - 0 • Email: [sales@medi-globe.de](mailto:sales@medi-globe.de)  
[www.medi-globe.de](http://www.medi-globe.de)

# SonoTip® ProControl EUS FNA – Single Use



## SonoTip® ProControl

With atraumatic (round) nitinol stylet; glass beads-blasted surface; incl. VacLok™ aspiration syringe (5–20 ml)

Catalog No.	Description	Needle Diameter (Gauge / mm)	Tube diameter (mm)	Needle length adjustment (cm)
GUS-33-18-025	EUS-Guided FNA System	25/0.5	1.8	0–8.5
GUS-33-27-025	EUS-Guided FNA System with tip stabilizer to accommodate large working channels (3.2 mm–3.8 mm)	25/0.5	2.7	0–8.5



## SonoTip® ProControl

With atraumatic (round) nitinol stylet; laser engraved surface; incl. VacLok™ aspiration syringe (5–20 ml)

Catalog No.	Description	Needle Diameter (Gauge / mm)	Tube diameter (mm)	Needle length adjustment (cm)
GUS-33-18-022	EUS-Guided FNA System compatible with guide wire max. .018"	22/0.7	1.8	0–8.5
GUS-33-27-022	EUS-Guided FNA System with tip stabilizer to accommodate large working channels (3.2 mm–3.8 mm) compatible with guide wire max. .018"	22/0.7	2.7	0–8.5

# SonoTip® ProControl EUS FNA – Single Use



## SonoTip® ProControl

With traumatic (beveled) nitinol stylet; laser engraved surface; incl. VacLok™ aspiration syringe (5–20 ml)

Catalog No.	Description	Needle Diameter (Gauge / mm)	Tube diameter (mm)	Needle length adjustment (cm)
GUS-34-18-022	EUS-Guided FNA System compatible with guide wire max. .018"	22/0.7	1.8	0–8.5
GUS-34-27-022	EUS-Guided FNA System with tip stabilizer to accommodate large working channels (3.2 mm–3.8 mm) compatible with guide wire max. .018"	22/0.7	2.7	0–8.5



## SonoTip® ProControl

With atraumatic (round) nitinol stylet; laser engraved surface; incl. VacLok™ aspiration syringe (5–20 ml)

Catalog No.	Description	Needle Diameter (Gauge / mm)	Tube diameter (mm)	Needle length adjustment (cm)
GUS-33-21-019	EUS-Guided FNA System compatible with guide wire max. .030" (recommended guide wire catalog no. GGW-08-30-400)	19/1.1	2.1	0–8.5