



Medi-Globe[®]
A MEDI-GLOBE COMPANY

ESD/EMR Knife

More Efficient and Economic



Used with an electrosurgical unit for endoscopic mucosal marking, incision in the GI tract.

ESD/EMR Knife

More Efficient and Economic

Catalog No.	Sales Unit	Outer Ø	Working Length	Requested min. Channel	Length I - Knife	Length O - Knife
MEK-401D-180	Box of 10	2.4 mm	180cm	2.8mm	2.0mm	4.0mm
MEK-402D-230	Box of 10	2.4 mm	230cm	2.8mm	2.0mm	4.0mm
MEK-403D-180	Box of 10	2.4 mm	180cm	2.8mm	3.0mm	4.0mm
MEK-404D-230	Box of 10	2.4 mm	230cm	2.8mm	3.0mm	4.0mm
MEK-405D-180	Box of 10	2.4 mm	180cm	2.8mm	4.0mm	4.0mm
MEK-406D-230	Box of 10	2.4 mm	230cm	2.8mm	4.0mm	4.0mm
MEK-407D-180	Box of 10	2.4 mm	180cm	2.8mm	1.0mm	4.0mm
MEK-408D-230	Box of 10	2.4 mm	230cm	2.8mm	1.0mm	4.0mm
MEK-409D-180	Box of 10	2.4 mm	180cm	2.8mm	1.5mm	4.0mm
MEK-410D-230	Box of 10	2.4 mm	230cm	2.8mm	1.5mm	4.0mm
MEK-411D-160	Box of 10	2.4 mm	160cm	2.8mm	2.0mm	4.0mm
MEK-412D-160	Box of 10	2.4 mm	160cm	2.8mm	3.0mm	4.0mm
MEK-413D-160	Box of 10	2.4 mm	160cm	2.8mm	4.0mm	4.0mm
MEK-414D-160	Box of 10	2.4 mm	160cm	2.8mm	1.0mm	4.0mm
MEK-415D-160	Box of 10	2.4 mm	160cm	2.8mm	1.5mm	4.0mm



- I knife–Designed for the purpose of marking, incision, peeling, hemostasis for mucosal/submucosal defects quickly.
- O knife–Ceramic tip designed to insulate for providing support to perform cutting and dissection more efficiency.
- 2 in 1 designed to accommodate physician technique preference, and more efficient, economic.

Medi-Globe GmbH

Medi-Globe-Straße 1-5
D-83101 Rohrdorf
OT Achenmühle

Tel: +49 (0) 8032 - 973-379
Fax: +49 (0) 8032 - 973-392
E-Mail: sales@medi-globe.de
www.medi-globe.de