

Spray Catheters

Single Use



**Application of Pigmented Liquids for Marking
of Suspicious Mucosa (Chromoendoscopy)**

Smooth Radiopaque Metal Tip

Atraumatic spray head

Spray Catheters – Single Use



Wide Spray



Radial Spray



Spray Catheter

Incl. stylet

Catalogue No.	Use for	Spray Head Diameter (mm)	Catheter Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
E47185-A	Wide Spray	1.8	1.8	120	2.0
E47185-C	Wide Spray	1.8	1.8	180	2.0
E47185-G	Wide Spray	1.8	1.8	230	2.0
E47226-C	Radial Spray	2.3	2.3	180	2.8
E47226-G	Radial Spray	2.3	2.3	230	2.8



LEF-GSC-SU-00-EN / LL-0360 Vers. 1.0 (08/2020)



ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
Email: sales@endo-flex.de • www.endo-flex.de