

Colonoscopy Accessories

Single Use



Suction Polyp Trap™

Aspiration containers for collection of polyps and specimens

Colonoscopy Shorts

Universal size with slit-on backside to protect patients' intimacy

Colon Decompression Probe

Radiopaque one-lumen catheter

Colonoscopy Accessoires – Single Use



Suction Polyp Trap™

For polyp collection after excision, sealable for transportation

Catalogue No.	Description
900334H	4 separated penetrable chambers to collect and store polyps in different sections, 4 open areas for straight unfiltered suction, 2 connection hoses at the lid



Colonoscopy Shorts

Consumables for instruments for medical purposes

Catalogue No.	Description
ZB0010	Three-ply material, water-repellent, bacteria-restraining, air permeable, universal size, slit-on backside



Colon Decompression Probe

Funnel connector, radiopaque tube, compatible with guide wire ≤ .052", non-sterile

Catalogue No.	Description	Outer Sheet Diameter (French)	Total Length (cm)
291140	Radiopaque Tube with 6 Distal Side Openings	14	175
291160	Radiopaque Tube with 6 Distal Side Openings	16	175
291180	Radiopaque Tube with 6 Distal Side Openings	18	175

LEF-GAC-SU-00-EN / LL-0374 Vers. 1.0 (10/2020)



ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
 Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
 Email: sales@endo-flex.de • www.endo-flex.de