

# Cytology Brushes

Single Use



## **Brush with Metal Tip**

Efficient tube closure to secure smear yield

## **Double Lumen Version**

Compatible with guide wire

# Cytology Brushes – Single Use



## Cytology Brush – Single Lumen

Catalogue No.	Working Length (cm)	Sheath Outer Diameter (mm)	Min. Accessory Channel (mm)
E4218-A	120	1.8	2.0
E4218-C	180	1.8	2.0
E4219-G	230	1.9	2.2
E4222-G	230	2.3	2.8



## Cytology Brush – Single Lumen

Catalogue No.	Working Length (cm)	Sheath Outer Diameter (mm)	Min. Accessory Channel (mm)
NE4218-A	120	1.8	2.0
NE4218-C	180	1.8	2.0
NE4222-G	230	2.3	2.8



## Cytology Brush – Double Lumen

Compatible with guide wire .035"

Catalogue No.	Working Length (cm)	Sheath Outer Diameter (mm)	Min. Accessory Channel (mm)
E4222DL-C	180	2.7	3.2



**ENDO-FLEX GmbH**

Alte Hünxer Str. 115 • 46562 Voerde • Germany  
 Tel.: +49 281 94400-0 • Fax: +49 281 94400-11  
 Email: sales@endo-flex.de • www.endo-flex.de