

# Pancreatic PTFE Stents

Single Use



## **Balanced Wall Diameters**

Maximum inner lumen with form stability and kink resistance

## **Tapered and Round Endings**

Atraumatic application


## **Multiple Holes**

Drainage of obstructed pancreatic ducts

## **Radiopaque Polytetrafluorethylen (PTFE) Material**

Controlled stent placement


# Pancreatic PTFE Stents – Single Use



## Pancreatic Stent – Curved, Perforated Throughout

Strong, incl. positioning sleeve; radiopaque; compatible with guide wire .035"

Catalogue No.	Size (French)	Length (cm)	Min. Accessory Channel (mm)
20505003	5	3	2.0
20505005	5	5	2.0
20505007	5	7	2.0
20505009	5	9	2.0
20507003	7	3	2.8
20507005	7	5	2.8
20507007	7	7	2.8
20507009	7	9	2.8



## Pusher

Compatible with guide wire .035"

Catalogue No.	Diameter Outer Tube (French)	Fitting to Stents (French)	Pusher Length (cm)
E20005001	5	5	180
E20007001	7	7	180



### ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany  
 Tel.: +49 281 94400-0 • Fax: +49 281 94400-11  
 Email: sales@endo-flex.de • www.endo-flex.de